

अखिलभारतीयआयुर्विज्ञानसंस्थान, रायपुर(छतीसगढ़)

All India Institute of Medical Sciences, Raipur (Chhattisgarh)

Tatibandh, GE Road, Raipur-492 099 (CG)

www.aiimsraipur.edu.in

Date:05/12/2019

AIIMS/R/2019/HS/MGPS/Oxygen Flow meter /LPC/155/

Inviting Quotations for Purchase of Oxygen flow meter with humidifier jar and BS type metal connector for the "Department of MGPS" AIIMS Raipur. For the MGPS Department AIIMS Raipur.

QUOTATION NOTICE

Sealed quotations are invited from intending registered Stockist / Distributors having TIN and relevant documents for Purchase of Oxygen flow meter with humidifier jar and BS type metal connector for the "Department of MGPS" AIIMS Raipur. For department of MGPS. The quotation with copy of certificate of TIN & other documents should be submitted to office of **Store Officer Room No. – 51, 2nd floor Ayush Buliding, Gate No. 1 up to 07/12/2019** before 12:00 pm. The quotations will be opened on the same day at 12:30pm. Details of item are given as under:-

Items Detail of Required Items									
S. no.	Complete Description of Items	Qty.	Make/ Model/ Brand	HSN Code	Basic rate price	GST %			ALL inclusive Cost
						SGST	CGST	IGST	
1.	Oxygen flow meter with humidifier jar and BS type metal connector	275							

Specification Back pressure compensated (BPC) flow meter with humidifier Jar:

Back Pressure Compensated (BPC) Flow meter made of brass body with or without humidifier with polycarbonate bottle is used for Oxygen therapy. The equipment Calibration Accuracy within \pm 5% of the set reading (International norm \pm 10% of the set reading) Inlet & Outlet nut is made leak proof by silicone O ring by hand tightening.

Flow indicating ball is made of Stainless Steel which spins steadily to ensure easy reading of flow.

1 low indicating ban is made of Stanness Steel which spins steading to ensure easy reading of now.			
Flow (L/Min)	0-15		
Increment (L/Min)	1		
Float	AISI 316 SS wire mesh		
Inlet & out let connection size	3/8" BSP RH		
Filter	Must be there		
Body	Brass (chrome plated)		
Calibrated pressure and Temperature Accuracy	50psig(3.45bar)at2 1degreeC/60psig(4.13bar)+10%		
	of the scale value or 0.2 L/M whichever is higher		

Humidifier valve settings	5-7psig
Inlet/out let nut nipple material	Brass

The Humidifier bottle can be fitted with the outlet of the BPC flow meter. The jar must be made of Polycarbonate and the cap is made of ABS plastic materials. Humidifier jar must be autoclaveble material.

Note: It should be compatible with gas out let of Beacon medaes and Penlon make gas out let and as well as with Mox regulator.

Terms & Condition

- 1. Firm to mention Manufactures name/Make/Brand name in their quotation.
- 2. Taxes, if any (Kindly mention in above table) should be clearly mentioned in the offer.
- 3. Document relating to registration of firm i.e. Tin number and relevant document should be submitted along with quotation.
- 4. Supply should be done within 10 days after Placement of PO.
- 5. Price should be FOR Destination basis (i.e. concerned department).
- 6. Payment will be released after certification from MGPS Department.
- 7. Quotation Name/No. and due date of opening must be mentioned on top of envelops.
- 8. LD @ 0.5% of delayed supply per week or part of week for delay of supply of material subject to maximum up to 10% of delayed supply to be deducted.
- 9. AIIMS Raipur reserves the right to place order for full or part quantity to one or more firms. The AIIMS, Raipur reserves the right to increase/decrease the number of required quantity.
- 10. All other terms & condition as per GCC applicable.
- 11. Material to be delivered at **Gate No-01**, **MGPS Department**,
- 12. Validity of quotation should be 90 days from the date of opening
- 13. Warranty: 1 year product replacement and technical repair service if any problem with product except physical damage.

V. Sitaramu Stores Officer (H) AIIMS, Raipur (C.G.)