



अखिलभारतीयआयुर्विज्ञानसंस्थान, रायपुर(छत्तीसगढ़)  
**All India Institute of Medical Sciences, Raipur (Chhattisgarh)**  
**Tatibandh, GE Road,**  
**Raipur-492 099 (CG)**  
[www.aiimsraipur.edu.in](http://www.aiimsraipur.edu.in)

AIIMS/R/HS/16994 /155/2020

Date: 05/03/2020

**Inviting Quotations for Purchase of consumable items (Blood Glucostrips) from Hospital Service for IPD/OPD at AIIMS Raipur.**

### **QUOTATION NOTICE**

Sealed quotations are invited from intending registered Stockist / Distributors having GST and relevant documents for Purchase of consumable items **(Blood Glucostrips) from Hospital Service for IPD/OPD at AIIMS Raipur**. The quotation with copy of certificate of GST & other documents should be submitted to office of **Store Officer Room No. – 51, 2nd floor AYUSH Building, Gate No. 1 up to 07/03/2020 before 12:00 pm**. The quotations will be opened on the same day at 012:30pm. Details of item are given as under:-

Sr. No	Specification	Quantity Required	Unit	HSN Code	Brand	Unit Rate	GST			Amount
							SG ST	IG ST	CGS T	
1	<b>Blood Glucostrips (50/pack)</b>	<b>1000</b>								

### **SPECIFICATION**

1. Accuracy- the industry as of the moment is allowed up to a plus minus 20% error.
2. Sample size-should not be more than 5 micro liters.\
3. Test time- should not be more than 10 seconds.
4. Interference- interference with commonly used durg and other agents along with effect on blood glucose measurement should be declared by manufacture.
5. Shelf life- shelf life of strip should be more than one year
6. Shelf life of strip after opening test strip should be individually packed or should be usable up to six month after first opening of strip box.
7. Coding- glucometer should be free from coding or calibration equipments in order to minimize errors at user level.
8. Pack size of strip should be available in minimum pack size of 50 strips.

9. Trigger mechanism- glucometer should start testing only when the required sample size has been absorbed by the strip in order to minimized wastage of strips.
10. Operating tempretuter should be 10-40 degree centigrade.
11. Measurement range- most meters read in a range between 0.6mmol/l and 33.3mmol/l (10mg/dl and 600mg/dl).
12. Battery- battery should be able to 1000 tests approximately and preferable should be AAA battery which is economical if not AAA then battery should be easily available and adorable.
13. Approvals- glucometer should be USFDA approval/ CE marked/ should have equivalent regulatory approval

Note; - the selected firm will have to provide glucometer (10% of quantity) free of charges for requirement of main hospital & centres etc.

### **TERM & CONDITION**

1. Firm to mention Make/Brand name in their quotation.
2. Taxes, if any (Kindly mention in above table) should be clearly mentioned in the offer.
3. Document relating to registration of firm i.e. GST and relevant document should be submitted along with quotation.
4. Products are certified from **CE / ISO** /as applicable, the Certificate to this effect should be attached.
5. Supply should be done within 15 days after Placement of PO.
6. Price should be FOR Destination basis (i.e. concerned department).
7. 100% Payment will be released after certification from Department of Pharmacy Store. **Quotation Name/No. and due date of opening must be mentioned on top of envelops.**
8. LD @ 0.5% of delayed supply per week or part of week for delay of supply of material subject to maximum up to 10% of delayed supply to be deducted.
9. AIIMS Raipur reserves the right to place order for full or part quantity to one or more firms. The AIIMS, Raipur reserves the right to increase/decrease the number of required quantity.
10. All other terms & condition as per GFR 2017.
11. Material to be delivered at **Pharmacy store, AIIMS Raipur.**
12. **Validity of the quotation should be 90 days from the date of opening.**
13. **Firm to submit documentary evidence in support of their claim for GST at the time of submission of bills.**
14. **Firm to submit sample, if required.**

**SO (H)**

**AIIMS, Raipur (C.G)**