



अखिलभारतीयआयुर्विज्ञानसंस्थान, रायपुर(छत्तीसगढ़)
All India Institute of Medical Sciences, Raipur (Chhattisgarh)
Tatibandh, GE Road,
Raipur-492 099 (CG)
www.aiimsraipur.edu.in

AIIMS/R/HS/URO/214/166/2020/PAC

Date: 28/10/2020

NOC

Sub:- Regarding procurement of Laproscopic Polymer Locking Clip Applier for Urology Department at AIIMS, Raipur under PAC.

The institute is in the process of purchase of Laproscopic Polymer Locking Clip Applier for Urology Department at AIIMS Raipur from M/s Sunshine Healthcare, Raipur. (OEM: M/s Teleflex Medical, USA) on proprietary basis. The proposal is submitted by Urology Department at AIIMS, Raipur and PAC Certifications are attached which is to be uploaded on website.

The above documents are being uploaded for open information to submit objection/comments, if any from any manufacturer regarding proprietary nature of the Item within 07 days of issued/uploading of the notification by reference No. AIIMS/R/HS/URO/214/166/2020/PAC. The comments should be sent to Store Officer (H), Room No. -146, C-1 Block, Lower Ground Floor, Gate No. 01, AIIMS, Raipur or through email at storesofficer.hp@aiimsraipur.edu.in on or before 04/11/2020 up to 3.00 PM failing which it will be presumed that any manufacturer/vendor having no comments to offer and case will be decided on merits.

Store Officer (H)
AIIMS Raipur (CG)

Encl:-

1. Proprietary letter of Vendor
2. Authentication letter of vendor
3. Certification for Purchase of Proprietary Article



अखिल भारतीय आयुर्विज्ञान संस्थान, रायपुर (छ.ग.)
All India Institute of Medical Sciences, Raipur (Chhattisgarh)
Tatibandh, GE Road,
Raipur-492 099 (CG)
www.aiimsraipur.edu.in

स्वामित्व प्रमाण पत्र
Proprietary Article Certificate

फाइल संख्या और संदर्भ File Number and Reference		No. 080109/2020/00/214 Date- 30-09-2020
1	सामाग्री का विवरण Description of article	Laparoscopic Polymer Locking Clip Applier (medium & Large size)
2	पूर्वानुमानित मात्रा/वार्षिक आवश्यकता Forecast of quantity/annual requirement	ONE SET
3	उपरोक्त मात्रा हेतु अनुमानित मूल्य Approximate estimated value for above	
4	निर्माता का नाम एवं पता Maker's name and address	TELEFLEX MEDICAL PVT. LTD. , NORTH CAROLINA, USA
5	अधिकृत डीलर/स्टाकिस्ट का नाम Name(s) of authorised dealers/stockists	SUNSHINE HEALTHCARE, RAIPUR (C.G.) - 492001
6	<p>मैं पी ए सी के आधार पर उपरोक्त खरीद को स्वीकार करता हूं और यह प्रमाणित करता हूं कि:</p> <p>नोट- (बी), (सी-1) या (सी-2) में से केवल एक को बनाए रखने के लिए टिक करें, जो भी लागू हो और दूसरो को काट दें। कृपया (ए) टिक कर पुष्टि करें इसके बिना पीएसी प्रमाण पत्र अवैध होगा</p> <p>I approve the above purchase on PAC basis and certify that:-</p> <p>Note- Tick to retain only one out of (b), (c-1) or (c-2) whichever is applicable and cross out others. Please do confirm (a) by ticking it – without which PAC certificate will be invalid.</p>	
6 (a)	यह एकमात्र फर्म है जो इस मद का निर्माण/संग्रहण कर रहा है। और This is the only firm who is manufacturing /stocking this item. AND	✓
6 (b)	किसी अन्य फर्म द्वारा समरूप मद निर्मित/विक्रय नहीं किया जाता है, जिसका उपयोग इसके बदले किया जा सकता है। अथवा A similar article is not manufacturing/sold by any other firm, which could be used in lieu OR	✓
6 (c-1)	कोई अन्य मेक/ब्रांड निम्नलिखित कारणों (जैसे ओईएम/वारंटी के) के लिए उपयुक्त नहीं होगा। अथवा No other make/brand will be suitable for following tangible reasons (like OEM/warranty spares): OR	—
6 (c)	कोई अन्य मेक/ब्रांड निम्नलिखित कारणों से उपयुक्त नहीं होगा (अगर पीएसी	—

	<p>पिछले खरीद में भी दिया गया था, तो कृपया इसके बाद से अधिक स्रोतों का पता लगाने के लिए प्रयास करें) तथा</p> <p>No other make/brand will be suitable for following intangible reasons (if PAC was also given in the last procurement cycle, please also bring out efforts made since then to locate more sources): OR</p>	—
	—
7	<p>प्रस्ताव के लिए वित्त शाखा की सहमति का संदर्भ (कार्रवाई भंडार और लेखा विभाग द्वारा की जायेगी)</p> <p>Reference of concurrence of finance wing to the proposal (Action will be taken by stores & Account Department)</p>

पिछले तीन सालों में इस मद की पीएसी खरीद का इतिहास नीचे दिया जा सकता है (यदि कोई हो) History of PAC purchase of this item for past three years may be given below (if any)

प्रदायक का नाम Name of the Supplier			
आदेश/निविदा संदर्भ और दिनांक Order/Tender reference & Date	आदेशित मात्रा Quantity Ordered	आदेश पर मूल दर (₹) Basic Rate on order (Rs.)	प्रतिकूल प्रदर्शन रिपोर्ट अगर कोई हो Adverse Performance Reported if any

अनुमोदन करने वाले प्राधिकारी का हस्ताक्षर -----

दिनांक 30-09-2020 -----

अधिकारी का पदनाम -----

Dr. Deepak Kumar Biswal

डॉ. दीपक कुमार बिस्वाल
Dr. Deepak Kumar Biswal
Assistant Professor (MS, MCh)
Department of Urology
All India Institute of Medical Sciences,
Raipur (C.G.) 492099

Technical Specification for Laparoscopic Polymer Locking-Clip Ligation System: (medium or Large size)

A. Clips should be made of inert, non-toxic, non-absorbable biocompatible polymer material. It should be radiolucent, nonconductive and does not interfere with CT, MRI or X-ray diagnostics.

B. There should be variety of ligating clip sizes for the ligation of vessels and/or tissues between 2 to 16 mm size.

C. The system should be compatible for use in laparoscopic surgeries.

D. Clips should be securely retained to applier jaws during application process.

E. Clips should be able to lock securely after applying to tissue/vessel for secure ligation.

F. Clip should have tactile feedback for confirmation of secure closure.

G. Clip ligating surface should have non-traumatic integrated teeth interface design to prevent slippage.

H. Clips should be arranged in an ergonomically designed cartridge for fast application.

I. If necessary, clips should allow removal with appropriate removal instrument.

J. Clip appliers should be made of corrosion free German steel with a matte finish to avoid glare during surgery.

K. Corresponding Laparoscopic Clip appliers should be colour-coded to match the colours of the corresponding Clip cartridges.

L. Laparoscopic Clip appliers should have 360° rotation knob facilitates for improved visibility & easy access to vessels.

M. All reusable Clip appliers should be compatible for sterilization with immersion, ETO & Sterrad method.

N. The system should be BIS/ European-CE (with four digit notified body number) & US-FDA approved product.

O. Should supply following Laparoscopic Polymer Locking-Clip Appliers:

Sr. No.	Clip Applier/ Remover	Shaft diameter	Working Length	Required Qty.
1.	Endo-surgery Clip Applier for M-size Clip	10mm	32cm	01
2.	Endo-surgery Clip Applier for L-size Clip	10mm	32cm	01

P. Warranty of Applier: One year.

Dr. Deepak Kumar Biswal
 Dr. Deepak Kumar Biswal
 Assistant Professor (MS, MCh)
 Department of Urology
 All India Institute of Medical Sciences,
 Raipur (C.G.) 492099

Teleflex®

Teleflex Medical
2917 Weck Drive
Research Triangle Park
North Carolina
27709
USA


Phone: +1-919-544-8000
Fax: 866-804-9881


Proprietary Certificate

To Whom it may Concern:

"Teleflex Medical", having office at 2917 Weck Drive, Research Triangle park, NC 27709 USA hereby confirms and verifies that following:

1. "Teleflex Medical" owns all the rights title and interest in the patents, trademarks, trade names, inventions, copyrights, know how, trade secrets and any other intellectual property of WECK HEM-O-LOK, the original Hem-o-lok invented and designed by Edward Weck, which falls under the Ligation Clip category for use in procedures involving ligation of vessels or tissue structures. Surgeon should apply appropriate size clip for the size of vessel or tissue structures to be ligated so that the clip completely encompasses the vessel or tissue structures. Weck Hem-o-lok ligating clips are available in a range of sizes starting from medium to extra-large for ligating vessel sizes from 2mm to 16mm. Hem-o-lok clip applicators are available in open surgery format as well as endoscopic surgery format with different sizes and angles based on surgeon's preference.
2. The following are the details of the /product:
 - Hem-o-lok ligating clips are made of acetal homopolymer that is inert, non-conductive, radiolucent and does not interfere with CT, MIR or X-ray diagnostics.
 - Hem-o-lok ligating clips have bosses designed to retain the clip in the applicator jaws.
 - Hem-o-lok ligating clips are bow-shaped allowing it to be removed with a Hem-o-lok clip remover, if needed.
 - Hem-o-lok ligating clips have an integrated teeth interface with the vessel and are designed to prevent slippage.
 - Hem-o-lok ligating clips have a locking mechanism that provides tactile feedback and secure closure.
 - Hem-o-lok clip applicators are reusable devices compatible with Weck Hem-o-lok ligating clips.
 - Hem-o-lok clip applicators are made from German steel with a matte finish to avoid any visual disturbance during surgery.
 - Hem-o-lok open clip applicators open at the box lock for easy, thorough cleaning.
 - Hem-o-lok endoscopic clip applicators have a 360-degree rotation knob for ease of use as well as a flush port for cleansing purposes.


Jennifer Perry
Sr. Regulatory Affairs Specialist


Dr. Deepak Kumar Biswal
Assistant Professor (MS, MCh)
Department of Urology
All India Institute of Medical Sciences,
Raipur (C.G.) 492099

Teleflex®

Teleflex Medical Private Limited

Registered Office ::

Building Embassy Vogue, 1st floor
Palace Road, Ward No. 78
Bangalore - 560 052, Karnataka, India
Phone: +91 80 40934790
Email: secretarial@teleflex.com
www.teleflex.com
CIN: U33112KA1993PTC124992

September 11, 2020

To
All India Institute of Medical Sciences,
Raipur

Dear Sir,

Subject: Authority letter for AIIMS Raipur - Supplies through Local Purchase


We, M/s Teleflex Medical Pvt Ltd, who are established and reputable manufacturers of Medical Disposables, having registered office at 1st Floor, Embassy Vogue Building, Palace Road, Bangalore - 560052, hereby authorize M/s. Sunshine Healthcare, Ward No. 51, Irrigation Colony, Near NHMM Hospital, Raipur, Chhattisgarh 492001 to collect orders, make supply, raise invoices and collect payment in their name on behalf of us.

We hereby extend our full guarantee and warranty as per the conditions for the goods offered for supply against this local purchase by the above firm.

The authorization is valid up to March 2021.

Yours faithfully,

For Teleflex Medical Pvt Ltd



Rajesh T.V

Director - Surgical & Govt Business

