

अखिल भारतीय आयुर्विज्ञान संस्थान ,रायपुर) छत्तीसगढ़ (All India Institute of Medical Sciences, Raipur (Chhattisgarh)

Tatibandh, GE Road, Raipur-492 099 (CG) <u>www.aiimsraipur.edu.in</u>

No. AIIMS/R/CS/BPS/189/21/LPC/A

Dated: 20/12/2021

Notice for inviting Quotation for supply of **Non Consumable items i.e. Dual Cuff Pneumatic Tourniquet for Department of Burn & Plastic Surgery, AIIMS, Raipur.**

OUOTATION NOTICE

Sealed quotations are invited from intending registered Stockiest/ Distributors having GST No./relevant documents for supply of Non Consumable item for Department of Burns & Plastic Surgery to be used in AIIMS, Raipur and should be submitted to Medical College Building, 2nd Floor, Gate no. 05 office of Store Officer up to 3:00 pm on 30.12.2021. The quotations will be opened on the same day at 3:30pm, item description is detailed bellow:

	Re	Make		HSN			
	~	Make	Unit	Code	Unit Price	GST @ %	Total Amount
ual Cuff Pneumatic Tourniquet	q. 01						
Dual cuffs: Dual port, color coded Dual display with integrated time coording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of ower allure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when tempt is made to deflate the machine Cuff lockout safety feature Cuff pressure: Adjustable 100 to 50 mmhg Pre-set resolution: 10mmhg Check Elapsed time: 0.000 - 3.000 ours Display: Three digit pre-selected pressure Three digit clock elapsed time Three digit clock elapsed time Toro Deflate Alarms: Cuff pressure is out of range For fore than 30sec.after initial start Cuff pressure is out of range for nore than 11 sec. during operation. After Elapsed time of 30min and 0min thereafter. See Three digits pre-selected pressure for the control of 30min and 0min thereafter.							
	Dual display with integrated time cording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of wer lure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when empt is made to deflate the machine Cuff lockout safety feature Cuff pressure: Adjustable 100 to 0 mmhg Pre-set resolution: 10mmhg Clock Elapsed time: 0.000 - 3.000 urs Display: Three digit pre-selected pressure is 3 digit clock elapsed time ii Green LED for inflate & red LED for Deflate. Alarms: Cuff pressure is out of range For one than 30sec.after initial start Cuff pressure is out of range for one than 11 sec. during operation. After Elapsed time of 30min and	Dual display with integrated time cording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of wer lure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when empt is made to deflate the machine Cuff lockout safety feature Cuff pressure: Adjustable 100 to 0 mmhg Pre-set resolution: 10mmhg Clock Elapsed time: 0.000 - 3.000 urs Display: Three digit pre-selected pressure i 3 digit clock elapsed time ii Green LED for inflate & red LED for Deflate. Alarms: Cuff pressure is out of range For one than 30sec.after initial start Cuff pressure is out of range for one than 11 sec. during operation. After Elapsed time of 30min and min thereafter. 3 beeps - 9 min.prior to set time	Dual display with integrated time cording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of wer lure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when empt is made to deflate the machine Cuff lockout safety feature Cuff pressure: Adjustable 100 to 0 mmhg Pre-set resolution: 10mmhg Clock Elapsed time: 0.000 - 3.000 urs Display: Three digit pre-selected pressure i 3 digit clock elapsed time ii Green LED for inflate & red LED for Deflate. Alarms: Cuff pressure is out of range For one than 30sec.after initial start Cuff pressure is out of range for one than 11 sec. during operation. After Elapsed time of 30min and min thereafter. 3 beeps - 9 min.prior to set time	Dual display with integrated time cording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of wer lure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when empt is made to deflate the machine .Cuff lockout safety feature .Cuff pressure: Adjustable 100 to 0 mmhg . Pre-set resolution: 10mmhg .Clock Elapsed time: 0.000 - 3.000 urs Display: .Three digit pre-selected pressure i 3 digit clock elapsed time ii Green LED for inflate & red LED for DeflateAlarms: Cuff pressure is out of range For one than 30sec.after initial start Cuff pressure is out of range for one than 11 sec. during operation. After Elapsed time of 30min and min thereafter. 3 beeps - 9 min.prior to set time	Dual display with integrated time cording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of wer lure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when empt is made to deflate the machine .Cuff lockout safety feature .Cuff pressure: Adjustable 100 to 0 mmhg . Pre-set resolution: 10mmhg .Clock Elapsed time: 0.000 - 3.000 urs Display: Three digit pre-selected pressure i 3 digit clock elapsed time ii Green LED for inflate & red LED for DeflateAlarms: Cuff pressure is out of range For ore than 30sec.after initial start Cuff pressure is out of range for ore than 11 sec. during operation. After Elapsed time of 30min and min thereafter. 3 beeps - 9 min.prior to set time	Dual display with integrated time cording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of wer lure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when empt is made to deflate the machine .Cuff lockout safety feature .Cuff pressure: Adjustable 100 to 0 mmhg . Pre-set resolution: 10mmhg .Clock Elapsed time: 0.000 - 3.000 urs Display: Three digit pre-selected pressure if 3 digit clock elapsed time if Green LED for inflate & red LED for DeflateAlarms: Cuff pressure is out of range For one than 30sec.after initial start Cuff pressure is out of range for one than 11 sec. during operation. After Elapsed time of 30min and min thereafter. 3 beeps - 9 min.prior to set time	Dual display with integrated time cording Use ambient air Microprocessor controlled Unit Self-check calibration Audible and Visual alarms in case of wer lure Internal pump for fast inflation time Positive locking connectors Alert the user of the cuff status when empt is made to deflate the machine . Cuff lockout safety feature . Cuff pressure: Adjustable 100 to 0 mmhg . Pre-set resolution: 10mmhg . Clock Elapsed time: 0.000 - 3.000 urs Display: Three digit pre-selected pressure i 3 digit clock elapsed time ii Green LED for inflate & red LED for Deflate . Alarms: Cuff pressure is out of range For one than 30sec.after initial start Cuff pressure is out of range for one than 11 sec. during operation. After Elapsed time of 30min and min thereafter. 3 beeps - 9 min.prior to set time

e) In deflate mode -cuff pressure is				
higher than 20. mmhg for 15 sec.after				
deflate mode is				
activated				
15. Dual compressor to allow for				
independent				
control over cuff one and cuff two				
16. Dedicated pressure line and				
pressure				
measurement line for cuff one and cuff				
two which is designed to provide most				
accurate				
readings				
17. Having Utilization of two				
independent cuffs during one surgical				
procedure to allow bilateral extremely				
procedures.				
18. Battery backup of two hours. So				
that it can be used during patient				
transport and to be used during power				
failure. No interruption in the				
procedure while the cuff is inflated -				
Automatically transfer power to the				
· · · · · · · · · · · · · · · · · · ·				
battery.				
19. Attached ability for 7 sizes of cuff				
from 8" to 42" of cylindrical cuffs				
meets individual requirements of thin				
and fat patient for arm and thigh cuffs.				
20. Cuffs will be of fabric with silicon				
bladder and bladder length equivalent				
to complete cuff length				
21. Cuffs will be autoclavable				
22. Will provide Double cuff for				
children				
23. Will provide Double cuff for long				
arm				
24. Will provide Double cuff for long				
1 1 -				
leg				
25. Will provide Double cuff for				
standard arm				
26. Will provide Double cuff for				
<u> </u>				
standard leg				
27. Instrument is BIS/ CE Certified.				
28. Company is ISO Certified.				
unit along with following cuffs,				
29. Warranty 1 year.				
29. Wairanty 1 year.		l		

Terms & Condition:

- 1. Rate should be mentioned in words & figure both.
- 2. Taxes/GST, if any (Kindly mention in table)
- 3. Delivery schedule: within 15 days from the date of issue of P.O.
- 4. No additional documents related to this NIQ will be entertained after opening of NIQ.
- 5. Price should be FOR Destination basis (i.e. Department of Burn & Plastic Surgery, AIIMS Raipur).

- 6. LD @ 0.5% of delayed supply per week or part thereof for delay of supply of material subject to maximum up to 10% should be deducted. After expiry of delivery period material cannot be accepted without extension of delivery period.
- 7. Quotation No/Name and due date of opening must be written on top of envelop.
- 8. Brand, Make & Warranty should be clearly mentioned in offer as well as tender/ quotation. Specific authorization may be submitted with the offer/bid.
- GST rates applicable on your quoted item may please be confirmed. HSN code for each item should be clearly mentioned.
- 10. Any change (Upward/Reduction) in the Basic Price structure post submission of bid shall be promptly intimated to AIIMS Raipur. You are also required to pass the Input Credit as per the following
- 11. Anti Profiteering Clause of GST.
 - "Upon Implementation of GST, any reduction in the rate of tax on supply of goods or service or the benefit of input tax credit shall be passed on to AHMS Raipur by way of commensurate reduction in the prices".
- 12. The GST registration details may be furnished.
- 13. RTGS details required for payment purpose.
- 14. 100% payment against receipt and acceptance of material.
- 15. No part supply or part payment will be entertained.
- 16. Validity of offer should not be less than 90 days.
- 17. Supply, installation and commissioning will be done by firm (if applicable).
- 18. The quantity in above column is totally tentative. It can be increased or decreased at the time of placement of order.
- 19. AIIMS, Raipur reserve the right to place the order for full or part quantity to one or more firms.
- 20. The product should be of superior quality and highly standard.

Store Officer, AIIMS, Raipur (C.G.)