

अखिलभारतीयआयुर्विज्ञानसंस्थान, रायपुर(छत्तीसगढ़) All India Institute of Medical Sciences, Raipur (Chhattisgarh) Tatibandh, GE Road, Raipur-492 099 (CG) www.aiimsraipur.edu.in

AIIMS/R/HS/SUR-GAS/119/2023/PAC

Call for Objections

Sub:- Inviting Comments/objections (if any), before declaring proprietary article for procurement of "Laparoscopic powered endocutter stapler and cartridge" for Surgical Gastroentrology Department at AIIMS Raipur.

Surgical Gastroentrology Department at AIIMS Raipur has to procure consumables for Laparoscopic powered endocutter stapler and cartridge through Proprietary Article basis under GFR -166.

The proposal submitted by M/s Johnson & Johnson Private limited India, Mumbai who is sole manufacturers and M/s Sai Associates, Durg is the local Agent of this item along with Proprietary Article Certificate are attached and same upload on AIIMS Raipur website.

The above documents are being uploaded for open information to submit objections/comments (if any) from any manufacturer/supplier before declaring proprietary article for the said equipment/items to be procured, within 07 days (i.e. 05/04/2023) from the date of issuance/uploading of the notifications.

The objection should be raised on the technical compliance sheet as enclosed, if any firm claiming suitability of their product with respect to specification mentioned.

The comments should be sent to the office of Stores Office (Hospital), Room No. 329, C C1 Block. No. 01. AIIMS. Raipur in sealed envelope through mail storesofficer.hp@aiimsraipur.edu.in before on 05/04/2023 up to 3.00 PM failing which it will be presumed that any other vendor having no comments to offer and case will be decided on merits.

Encl:-

- 1. Proprietary Letter by Manufacturer
- 2. Proprietary Letter by Department
- 3. Authentication letter of vendor

Stores Officer- H AIIMS Raipur

Date: 29/03/2023

भंडार आधकारी (अस्प.) Store Officer (H) एम्स, रायपुर (छ.ग.) AIIMS, Ralpur (C.C)



अखिलभारतीय आयुर्विज्ञानसंस्थान, रायपुर (छ.ग.) All India Institute of Medical Sciences, Raipur (Chhattisgarh) Tatibandh, GE Road, Raipur-492 099 (CG)

Raipur-492 099 (CG) www.aiimsraipur.edu.in

स्वामित्वप्रमाणपत्र Proprietary Article Certificate

| | संख्या औरसंदर्भ | | |
|--------|--|--|--|
| ile Nu | mber and Reference | | |
| 1 | सामाग्रीकाविवरण | Laparoscopic powered endocut stapler and cartridge | |
| | Description of article | stapler and carriage | |
| 2 | पूर्वानुमानितमात्रा / वार्षिकआवश्यकता | | |
| | Forecast of quantity/annual requirement . | | |
| 3 | उपरोक्तमात्रा हेतुअनुमानितमूल्य | | |
| | Approximate estimated value for above | | |
| 4 | निर्माताकानाम एवं पता | | |
| | | M/s Johnson & Johnson Pvt Ltd, India | |
| 5 | अधिकृतडीलर / स्टाकिस्टकानाम | M/S Sai Associates, Durg | |
| 3 | Name(s) of authorised dealers/stockists | | |
| | मैंपी ए सी के आधारपरउपरोक्त खरीदकोस्वीकारकरताहूं और यह | इ प्रमाणित करताहूंकिः | |
| | नोट- (बी), (सी-1) या (सी-2) में से केवल एक कोबनाए | रखने के लिए टिककरें, | |
| | जोभीलागूहोऔरदूसरोकोकाटदें।कृपया (ए) टिककरपुष्टिकरेंइसवे | विनापीएसीप्रमाणपत्र अवैध | |
| 6 | होगा | | |
| | I approve the above purchase on PAC basis and certify that:- | | |
| | Note- Tick to retain only one out of (b), (c-1) or (c-2) whichever is applicable and | | |
| | Note- Tick to retain only one out of (b), (c-1) or (c-2) whi | ichever is applicable and | |
| | Note- Tick to retain only one out of (b), (c-1) or (c-2) which cross out others. Please do confirm (a) by ticking it – with | | |
| | cross out others. Please do confirm (a) by ticking it – without will be invalid. | | |
| | cross out others. Please do confirm (a) by ticking it - with | | |
| 6 (a) | cross out others. Please do confirm (a) by ticking it – without will be invalid. | | |
| 6 (a) | cross out others. Please do confirm (a) by ticking it – withowill be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण/संग्रहणकररहाहै। | out which PAC certificate | |
| 6 (a) | cross out others. Please do confirm (a) by ticking it – withould will be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै। और This is the only firm who is manufacturing /stocking this item AND | n. | |
| 6 (a) | cross out others. Please do confirm (a) by ticking it – withould be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै। और This is the only firm who is manufacturing /stocking this iter. | n. | |
| | cross out others. Please do confirm (a) by ticking it – withould will be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै। और This is the only firm who is manufacturing /stocking this item AND | n. | |
| 6 (a) | cross out others. Please do confirm (a) by ticking it – withould will be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै। और This is the only firm who is manufacturing /stocking this iter AND किसीअन्य फर्म द्वारासमरूपमदनिर्मित / विक्रय नहीं कियाजाताहै, | n. | |
| | cross out others. Please do confirm (a) by ticking it – withowill be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै। और This is the only firm who is manufacturing /stocking this iter AND किसीअन्य फर्म द्वारासमरूपमदनिर्मित / विक्रय नही कियाजाताहै, जिसकाउपयोगइसकेबदलेकियाजासकताहै।अथवा | n. | |
| | cross out others. Please do confirm (a) by ticking it – withould will be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै। और This is the only firm who is manufacturing /stocking this iter AND किसीअन्य फर्म द्वारासमरूपमदनिर्मित / विक्रय नही कियाजाताहै, जिसकाउपयोगइसकेबदलेकियाजासकताहै।अथवा A similar article in not manufacturing/sold by any other firm | n. | |
| 6 (b) | cross out others. Please do confirm (a) by ticking it – withould will be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै । और This is the only firm who is manufacturing /stocking this item AND किसीअन्य फर्म द्वारासमरूपमदनिर्मित / विक्रय नही कियाजाताहै, जिसकाउपयोगइसकेबदलेकियाजासकताहै ।अथवा A similar article in not manufacturing/sold by any other firm which could be used in lieu OR | n. | |
| | cross out others. Please do confirm (a) by ticking it – with will be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै । और This is the only firm who is manufacturing /stocking this iter AND किसीअन्य फर्म द्वारासमरूपमदनिर्मित / विक्रय नही कियाजाताहै, जिसकाउपयोगइसकेबदलेकियाजासकताहै ।अथवा A similar article in not manufacturing/sold by any other firm which could be used in lieu OR कोईअन्य मेक / ब्रांडनिम्नलिखितकारणो (जैसेओईएम / वारंटी व | n. | |
| 6 (b) | cross out others. Please do confirm (a) by ticking it – withould will be invalid. यह एकमात्र फर्महैजोइसमदकानिर्माण / संग्रहणकररहाहै। और This is the only firm who is manufacturing /stocking this iter AND किसीअन्य फर्म द्वारासमरूपमदनिर्मित / विक्रय नहीं कियाजाताहै, जिसकाउपयोगइसकेबदलेकियाजासकताहै।अथवा A similar article in not manufacturing/sold by any other firm which could be used in lieu OR कोईअन्य मेक / ब्रांडनिम्नलिखितकारणों (जैसेओईएम / वारंटी व लिए उपयुक्त नहीं होगा।अथवा | n. | |

Receipt No. 9897/2023/SG (Computer No. 9897)

| 6 (c) | कोईअन्य मेक / ब्रांडनिम्नलिखितकारणों से उपयुक्त नहीं होगा (अगर पीएसीपिछले खरीदमेंभीदियागयाथा, तोकृपयाइसकेबाद से अधिकस्रोतोका पता लगाने के लिए प्रयास करें)तथा No other make/brand will be suitable for following intangible reasons (if PAC was also given in the last procurement cycle, please also bring out efforts made since then to locate more sources): OR |
|-------|--|
| | प्रस्ताव के लिए वित्त शाखा की सहमतिकासंदर्भ(कार्रवाई भंडारऔरलेखाविभाग द्वारा की जायेगी) Reference of concurrence of finance wing to the proposal (Action will be taken by stores & Account Department) |

| पिछलेतीनसालोमेंइसमद की पीए PAC purchase of this item for past | सी खरीदकाइति three years may | ाहासनीचेदियाजास v be given below (if | कताहै (यदि कोई हो)History of any) |
|--|-------------------------------------|---|---|
| प्रदायककानाम Name of the Supplier आदेश / निविदासंदर्भऔरदिनांक Order/Tender reference & Date | आदेशितमात्रा Quantity Ordered | आादेशपरमूलदर (□) Basic Rate on order (②) | प्रतिकूलप्रदर्शनरिपोर्टअगरकोईहो Adverse Performance Reported if any |
| none | | | |

Technical Specification

Laparoscopic powered endocutter stapler and cartridge - 1 set

It consists of:

| serial no | item | description | Quantity needed |
|--------------|--|--|--------------------|
| 1 | Battery Operated 60mm Articulating Endo Cutter | It should be Battery Operated 60mm Articulating Endo Cutter with a disposable battery pack; should have enhanced distal tip stability while firing; should have closed channel in the cartridge jaw for better stability during firing; 360 degree rotation shaft and one handed natural articulation up to 45 degrees; should have Precision Machined Anvil to deliver initial, system wide compression; should have Wide Proximal to Distal Jaw Aperture (proximal 8mm, Distal 22mm); should have 3 Point Gap Control for alignment and calibration throughout the 60 mm staple line; should have Knife direction/reverse control to discontinue the firing and return the knife; should allow interchangeable 6 row cartridge options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm and all fits down to 12mm trocar sleeve; 340 mm shaft length. (BIS/USFDA/CE certified as applicable) | 1 |
| 2 | Battery Operated 45mm Articulating Endo Cutter | It should be Battery Operated 45mm Articulating Endo Cutter with a disposable battery pack; should have enhanced distal tip stability while firing; should have closed channel in the cartridge jaw for better stability during firing; 360 degree rotation shaft and one handed natural articulation up to 45 degrees; should have Precision Machined Anvil to deliver initial, system wide compression; should have Wide Proximal to Distal Jaw Aperture (proximal 8mm. Distal 22mm); should have 3 Point Gap Control for alignment and calibration throughout the 60 mm staple line; should have Knife direction/reverse control to discontinue the firing and return the knife; should allow interchangeable 6 row cartridge options of closed staple height of 1.0, 1.5, 1.8, 2.0 and 2.3 mm and all fits down to 12mm trocar sleeve; 340 mm shaft length. (BIS/USFDA/CE certified as applicable) | 1 |
| 3 | Endoscopic linear cutter cartridge with pocket extension having 60mm staple line with 6 rows of staples, for regular thick tissues - 60 mm (gold) Endoscopic linear cutter cartridge with pocket extension having 60mm staple line with 6 rows of staples, for regular thick tissues, cartridge having closed staple height of 1.8 mm compatible with the above mentioned Battery operated 60 mm articulating endo cutter. (BIS/USFDA/CE certified as applicable) | | 5 |
| 4 | Endoscopic linear cutter cartridge for thick tissues - 60 mm (green) | Endoscopic linear cutter cartridge with pocket extension having 60mm staple line with 6 rows of staples, for thick tissues, cartridge having closed staple height of 2.0 mm compatible with the above mentioned Battery operated 60 mm articulating endo cutter. (BIS/USFDA/CE certified as applicable) | 5 |

Terms and conditions:

- 1. All the items should be from the same manufacturing company
- 2. Item wise selection as L-1 shall not be done
- 3. The bidder who has not quoted for all the instruments i.e. the entire set shall not be entertained and bid shall be rejected
- 4. Items must be USFDA, CE, ISO certified.

#Note 1: final quantity of cartridges may vary depending on the final rate by the bidding firm.



ETHICON PARTO PAR Columna Columna FARITY OF COMPANIES

Ethicon Endo Surgery LEG 475 Calle C Guaynabo, Puerto Rico 00969 USA

PROPRIETARY CERTIFICATE & AGREEMENT FOR REPRESENTATION

We hereby certify that product attached in Annexure A, are proprietary product of **Ethicon** Endo-Surgery, LLC part of Johnson & Johnson Family of Companies. Ethicon Endo-Surgery, LLC is the legal manufacturer of the products mentioned under Annexure A.

We further certify that M/s Johnson & Johnson Private Limited, India having registered office at L.B.S Marg, Mulund (W), Mumbai-400080 is our affiliate and is authorized to market the product in India through its regional branch office at 501, Arena Space, Jogeshwari (E) - 400060, India. Johnson & Johnson Private Limited, India and its regional branch office is authorized to submit offers and commit, sign, finalize, execute orders and also to appoint distributors and dealers for selling the complete range of products on its own behalf.

We also certify of the adequacy of technical expertise of the products offered by Johnson & Johnson private limited India and will give them back-up support toward supply, spares & technical updated.

Signed for and on behalf of: Ethicon Endo-Surgery, LLC

Brian Godwin

Associate Director - Regulatory Affairs

February 17, 2022



2 FEB 2023

ETHICON

PARTOR THE Johnson Johnson FAMILY OF COMPANIES

Ethicon Endo-Surgery, LLC 475 Caile C Guaynabo, Puerto Rico (XI989 USA

| ATS45NK | ENDOPATH ETS-Flex45 No Knife Endoscopic Articulating Linear Stapler | | |
|----------|---|--|--|
| SC60A | ECHELON FLEX 60 ENDOPATH Stapler Compact Articulation Endoscopic Linear Cutter | | |
| PLE60A | ECHELON FLEX 60 Powered ENDOPATH Stapler Long Articulating Endoscopic Linear Cutter | | |
| PSE60A | ECHELON FLEX 60 Powered ENDOPATH Stapler Articulating Endoscopic Linear Cutter | | |
| 2B12LT | ENDOPATH XCEL WITH OPTIVIEW TECHNOLOGY Bladeless Trocar with Stability Sleeve, 12mm, 100mm shaft length | | |
| 2B12XT | ENDOPATH XCEL WITH OPTIVIEW TECHNOLOGY Bladeless Trocar with Stability Sleeve, 12mm, 150mm shaft length | | |
| 2CB12LT | ENDOPATH XCEL WITH OPTIVIEW TECHNOLOGY Universal Trocar Stability Sleeve, 12mm, 100mm shaft length | | |
| ECS33A | Ethicon Endo-Surgery Endoscopic Curved Intraluminal Staplers (ILS) (33mm) | | |
| PCE45A | ECHELON FLEX 45 Powered ENDOPATH Stapler Compact Articulating Endoscopic Linear Cutter | | |
| PVE35A | ECHELON FLEX Powered Vascular Stapler with Advanced Placement Tip | | |
| VASECR35 | ENDOPATH ECHELON Vascular Reload | | |
| GST60W | ENDOPATHECHELON Endoscopic Linear Cutter Reload with Gripping Surface Technology, white – 2,6 mm | | |
| GST60B | ENDOPATHECHELON Endoscopic Linear Cutter Reload with Gripping Surface Technology, blue – 3,6 mm | | |
| GST60D | ENDOPATHECHELON Endoscopic Linear Cutter Reload with Gripping Surface Technology, gold – 3,8 mm | | |
| GST60G | ENDOPATH ECHELON Endoscopic Linear Cutter Reload with Gripping Surface Technology, green – 4,1 mm | | |
| GST60T | ENDOPATH ECHELON Endoscopic Linear Cutter Reload with Gripping Surface Technology, black – 4,2mm | | |
| PLEE60A | ECHELON FLEX 60 Powered Plus ENDOPATH Stapler Long Articulating Endoscopic Linear Cutter | | |
| PSEE60A | ECHELON FLEX 60 Powered Plus ENDOPATH Stapler Articulating Endoscopic Linear Cutter | | |
| EPH02 | ENDOPATH Electrosurgery Probe Plus II Pistol Grip Handles, Hand Control | | |
| EPS01 | ENDOPATH Electrosurgery Probe Plus II, 5 mm Hook, 34 cm length | | |
| EPS02 | ENDOPATH Electrosurgery Probe Plus II, 5 mm Spatula, 34 cm length | | |
| EPS03 | ENDOPATH Electrosurgery Probe Plus II, 5 mm Right Angle, 34 cm length | | |
| EPS04 | ENDOPATH Electrosurgery Probe Plus II, 5 mm Curved Dissector, 34 cm length | | |
| EPS05 | ENDOPATH Electrosurgery Probe Plus II, 5 mm Hook, 29 cm length | | |
| EPS06 | ENDOPATH Electrosurgery Probe Plus II, 5 mm Spatula, 29 cm length | | |
| EPS07 | ENDOPATH Electrosurgery Probe Plus II, 5 mm Right Angle, 29 cm length | | |
| GEN11 | Generator G11 | | |
| CRT11 | Generator G11 Cart | | |
| FSW11 | Footswitch G11 | | |
| HDH05 | HARMONIC Dissecting Hook, 5 mm diameter, 32cm length | | |
| HP054 | HARMONIC Hand Piece | | |
| HPBLUE | HARMONIC Blue Hand Piece | | |
| HAR23 | Harmonic ACE 5 mm Diameter Shears 23 cm Length + Adaptive Tissue Technology | | |

Receipt No. 9897/2023/SG (Computer No. 9897)

ETHICON

Ethicon Endo Surgery, LLC 475 Calle C Guaynaho, Puerto Rico 00959 USA

| PSEE45A | ECHELON FLEX 45 Powered Plus ENDOPATH Stapler Articulating Endoscopic Linear |
|----------|--|
| | Cutter |
| PCEE60A | ECHELON FLEX 60 Powered Plus ENDOPATH Stapler Compact Articulating Endoscopic Linear Cutter |
| GST45B | ENDOPATH ECHELON Endoscopic Linear Cutter Reload, 45mm, with Gripping Surface Technology, blue – 3,6 mm |
| GST45D | ENDOPATHECHELON Endoscopic Linear Cutter Reload, 45mm, with Gripping Surface Technology, gold – 3,8 mm |
| GST45G | ENDOPATH ECHELON Endoscopic Linear Cutter Reload, 45mm, with Gripping Surface Technology, green – 4,1 mm |
| GST45T | ENDOPATH ECHELON Endoscopic Linear Cutter Reload, 45mm, with Gripping Surface Technology, black – 4,2 mm |
| GST45W | ENDOPATH ECHELON Endoscopic Linear Cutter Reload, 45mm, with Gripping Surface Technology, white – 2,6 mm |
| HARHD20 | HARMONIC HD 1000i Shears, 5mm diameter, 20cm shaft length |
| EPH04 | ENDOPATH Electrosurgery Probe Plus II Pencil Grip Handles, Hand Control |
| NSLX120L | ENSEAL X1 Large Jaw Tissue Sealer |

Johnson Johnson

MEDICAL DEVICES

DIVISION OF JOHNSON & JOHNSON PRIVATE LIMITED

Dated: 28th Jan, 2023

To, The Director, All India Institute of Medical Sciences, Raipur (C.G)

Subject: Authority Letter to supply Ethicon range of the products.

Dear Sir,

This is to certify that M/s Sai Associates is out authorized distributor for Ethicon (A company of Johnson and Johnson) range of products for AIIMS Raipur for quoting prices, receiving orders, supply and collect payment on our behalf.

The contact details of M/S. Sai Associates are Sai Associates
Regd Off – E – 12, Om Complex, Arya Nagar,
Durg – 491001
CG (India)
Contact No – 7028180201

Assuring you of our best services always.

Yours truly

Saurabh Pandey Digitally signed by Saurabh Pandey DNI: cn=Saurabh Pandey, o, ou, email=spandey6eits.jnj.com, c=US Reason: I am the author of this document. Date: 2023.01.28 10:20:54 +05'30'

Saurabh Pandey

Regional Sales Manager – Ethicon Johnson & Johnson Medical india



Johnson House, 64/66, Senapati Bapat Marg, Post Box No. 16480, Mahim, Mumbai 400 016 India T: +91 22 24466801 / 24466818 F: +91 22 24449105 / 24448031

Registered Office: Johnson & Johnson Private Limited, L.B.S. Marg, Mulund (W), Mumbai 400 080.

T: +91 -22-2564 4441 F: +91-22-2564 3410

Web: www.jnjindia.com CIN: U33110MH1957PTC010928