

अखिल भारतीय आयुर्विज्ञान संस्थान, रायपुर(छत्तीसगढ़)

All India Institute of Medical Sciences, Raipur (Chhattisgarh) Tatibandh, GE Road, Raipur-492 099 (CG)

www.aiimsraipur.edu.in

Date: 25-04-2020

AIIMS/R/HS/G/2020/155/BMoT/58

Inviting Quotations for procurement of Blood and Fluid Warmer for Super Speciality B-Block MoT at AIIMS, Raipur.

QUOTATION NOTICE

Sealed quotations are invited from intending registered Stockist / Distributors having GST and relevant documents for procurement of Blood and Fluid Warmer for Super Speciality B-Block MoT at AIIMS, Raipur. The quotation with copy of certificate of GST & other documents should be submitted to office of Hospital Stores, Room No. – 146, C-1 Block, Opp. Department of Nuclear Medicine, Gate No.01, AIIMS, Raipur or through email on storesofficer.hp(@aiimsraipur.edu.in up to 30-04-2020 before 03:00 pm. The quotations will be opened on the same day at 03:30 pm. Details of item are given as under:-

Sl.	Name of Item	Quantity	Unit	HSN	Brand/	Unit	G	ST	Amount
01	Blood and Fluid Warmer Specification – As per Annexure-I	01	Nos						
	Grand Total				I		<u> </u>	<u> </u>	

Terms & Condition:

- 1. Firm to mention Make/Brand name in their quotation.
- 2. Taxes, if any (Kindly mention in above table) should be clearly mentioned in the offer.
- 3. Document relating to registration of firm i.e. GST and relevant document should be submitted along with quotation.
- 4. Products are certified from ISI/ISO/CE/BIS/CE/USFDA as applicable, the Certificate to this effect should be attached.
- 5. Supply should be done within 15 days after Placement of PO.
- 6. Price should be FOR Destination basis (i.e. concerned department).
- 7. 100% Payment will be released after certification from concerned department.
- 8. Quotation Name/No. and due date of opening must be mentioned on top of envelops.
- 9. LD @ 0.5% of delayed supply per week or part of week for delay of supply of material subject to maximum up to 10% of delayed supply to be deducted.
- 10. AIIMS Raipur reserves the right to place order for full or part quantity to one or more firms. The AIIMS, Raipur reserves the right to increase/decrease the number of required quantity.
- 11. All other terms & condition as per GCC applicable.

- 12. Material to be delivered at Super Speciality B-Block MoT 4th Floor, Raipur.
- 13. Validity of the quotation should be 90 days from the date of opening.
- 14. Firm to submit documentary evidence in support of their claim for GST at the time of submission of bills.
- 15. Firm to submit sample if required.

V. Sitaramu Stores Officer-H AIIMS, Raipur (CG)

Blood & Fluid Warmer - Smiths Medical / Hotline/ HL 90

Technical Specification-

Standard Compliance	Guidelines						
Product safety	EN60601-1, UL 2601-1						
EMC	EN 60601-1-2, FCC 47 CFR Part 15, Class B						
Enclosure Protection	IEC60529 IP Code : IPX1						
Fluid Warmers	ASTM F2172-02						
Physical							
Height, Overall	24.1 cm (9.5 inches)						
Width, Overall	21.0 cm (8.3 inches)						
Depth, Overall	17.8 cm (7.0 inches)						
Weight, Dry	3.5 Kg (7.6 lbs)						
Weight, Wet (with recirculating solution)	5.0 Kg (11.0 lbs)						
Weight, Shipping	3.6 Kg (7.95 lbs)						
Recirculating solution Capacity	1.4 L (0.37 gallons)						
Maximum Height on I.V. Pole	107 cm (42 inches)						
	207 cm (12 menes)						
Environmental	Temperature	Unami dia					
Operation	10°C to 45°C	Humidity					
Transportation	-18°C to 60°C	10 to 95					
Storage	-18°C to 60°C	5 to 90					
		5 to 90					
Thermal	Temperature						
Temperature Set Point	41.9_{0} C ± 0.1_{0} C						
Over Temperature Set Point	43.1 ₀ C						
Floatwicel	S ESCENIENCES						
Electrical MAINS Page 1	Туре						
MAINS Power Input 230V							
	230VAC, 50/60 HZ,1.5 Amps						
MAINS Auxiliary Supply Power Output:							
230V	220140 50 60 15						
	230VAC, 50/60 HZ,0.6 Amps						
Electrical	Tyme						
Protection Against Electrical Shock	Type						
Model Of Operation	Class 1 Equipment, Type BF						
Type of Current	Continuous						
ngress Protection Rating	Alternating IPX1						
Mulig	ILVI						
Performance							
Recirculating Solution Temperature	Recirculating solution temperature reaches 37°C from ambient in about 4 minutes						
Flow Rates	At gravity flow rates to 5,000 ml per hour						
Varrenty	01 Year						