



अखिल भारतीय आयुर्विज्ञान संस्थान, रायपुर (छत्तीसगढ़)  
All India Institute of Medical Sciences, Raipur (Chhattisgarh)  
Tatibandh, Raipur-  
492 099 (CG)

AIIMS /Micro/2025-26/106/ LPC/A/E-20936/1008

Date: 04/10/2025

विषय/Sub: Inviting Quotation for two bid system (Technical and Financial) for the procurement of "Consumables Items (1st Call)" at Microbiology Department of AIIMS Raipur.

## कोटेशन सूचना QUOTATION NOTICE

Two separate sealed quotation is invited from intending registered stock list/ Distributors having GST Number and relevant documents for "Consumables Items (1st Call)" for Microbiology Department, AIIMS Raipur. The quotation should be submitted to Department of Microbiology office, room no. 1220, Medical college building, first floor, Gate No. 5, AIIMS Raipur up to 03:00 pm on 10/10/2025. The quotations will be opened on the same day at 3:30 pm. Details of item are given as under:

| Sr. No. | NAME                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | Require ment | Unit | Make & Hsn Code | Unit Rate In Rs. | Gst In % | Unit Rate With Gst | Total Price In Rs. |
|---------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|------|-----------------|------------------|----------|--------------------|--------------------|
| 1.      | Multiplex PCR Panel for Gastro Intestinal Pathogens Gastro Pathogen Detection Kit Pack Size: 24 Reactions BACTERIA<br>• <i>Clostridium difficile</i> toxin A/B (tcdA/tcdB)<br>• Enterotoxigenic <i>Escherichia coli</i> (EAEC) (aatA/aggR)<br>• Enteroinvasive <i>Escherichia coli</i> (EIEC)/ <i>Shigella</i> (ipaH)<br>• Enteropathogenic <i>Escherichia coli</i> (EPEC) (eaeA)<br>• Enterotoxigenic <i>Escherichia coli</i> (ETEC) (eltA/estA/estB)<br>• <i>Helicobacter pylori</i><br>• Pathogenic <i>Campylobacter</i> spp ( <i>C. jejuni</i> , <i>C. upsaliensis</i> , <i>C. coli</i> )<br>• <i>Plesiomonas shigelloides</i><br>• <i>Salmonella</i> spp<br>• Shiga toxin-producing <i>E. coli</i> (STEC) (stx1/stx2)<br>• Shiga toxin-producing <i>E. coli</i> (STEC) O157:H7<br>• <i>Vibrio cholerae</i><br>• <i>Vibrio parahaemolyticus</i> | 01           | Kit  |                 |                  |          |                    |                    |

- *Vibrio vulnificus*
  - *Yersinia Enterocolitica*
- VIRUSES**
- Adenovirus (F40/41)
  - Astrovirus
  - Norovirus genogroup I&II
  - Rotavirus A
  - Sapovirus (I, II, IV, V)
- PROTOZOA**
- *Cyclospora cayetanensis*
  - *Cryptosporidium* spp
  - *Entamoeba Histolytica*
  - *Giardia Lamblia*

2. **Respiratory Bacterial**  
**Panel 50Rx**  
Pack Size: 50 Reactions  
**PATHOGENS COVERED**
- *Bordetella pertussis*
  - *Staphylococcus aureus*
  - *Klebsiella pneumoniae*
  - *Chlamydophila pneumoniae*
  - *Haemophilus influenzae*
  - *Haemophilus influenza B*
  - *Mycoplasma pneumoniae*
  - *Corynebacterium diphtheriae*
  - *Streptococcus pyogenes* (Group-A)
  - *Streptococcus pneumoniae*
  - *Group C and G streptococci*
  - *Corynebacterium ulcerans*
  - *Exogenous Control*
  - *Cryptococcus neoformans*
  - *Pneumocystis jirovecii*
  - *Aspergillus* species

3. **Multiplex Panel for Fungal Respiratory Pathogens**  
**Fungal Respiratory Panel TaqMan PCR Kit**  
Pack Size: 50 Reactions  
**PATHOGENS COVERED**  
Panfungal,  
Aspergillus,  
Mucorales,  
*Pneumocystis jirovecii*

01

Kit

01

Kit


|     |                                                                                                                                                                                                                                                                       |    |     |  |  |  |  |  |
|-----|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----|-----|--|--|--|--|--|
| 4.  | <u>Multiplex Panel for Tropical Fever Pathogens</u><br>Tropical Fever Panel (8-Pathogen) Panel TaqMan PCR Kit<br>Pack Size: 50 Reactions<br><b>PATHOGENS COVERED</b><br>Dengue<br>Chikungunya<br>Typhoid<br>Malaria<br>Leptospirosis<br>Scrub typhus,<br>Zika,<br>JEV | 01 | Kit |  |  |  |  |  |
| 5.  | Nucleic Acid Isolation Kit (Body Fluids)<br>Pack Size: 50 Reactions                                                                                                                                                                                                   | 02 | Kit |  |  |  |  |  |
| 6.  | Ultra Clean Tube Conical Bottom, with HDPE Closure, 15 sterile, Pack Size: 25x20                                                                                                                                                                                      | 1  |     |  |  |  |  |  |
| 7.  | Ultra Clean Tube Conical Bottom, with HDPE Closure, 50 sterile, Pack Size: 25x20                                                                                                                                                                                      | 1  |     |  |  |  |  |  |
| 8.  | Self-Standing Centrifuge Tube, with HDPE Closure<br>Pack Size: 25x20                                                                                                                                                                                                  | 1  |     |  |  |  |  |  |
| 9.  | Universal Combi Rack, Pack Size: 4 No.                                                                                                                                                                                                                                | 2  |     |  |  |  |  |  |
| 10. | Combilock Rack, Pack Size: 4 No.                                                                                                                                                                                                                                      | 2  |     |  |  |  |  |  |
| 11. | Vial External Threaded, Star-Foot, Sterile, 1.8 ml, Pack Size: 1000 No. per pack                                                                                                                                                                                      | 1  |     |  |  |  |  |  |
| 12. | Cryo Vial Workstation Rack - 50 Places, PC, Pack Size: 1                                                                                                                                                                                                              | 1  |     |  |  |  |  |  |
| 13. | Cryo Cube Box, Pack Size: 4                                                                                                                                                                                                                                           | 2  |     |  |  |  |  |  |
| 14. | Cryo Cube Box, Pack Size: 8                                                                                                                                                                                                                                           | 1  |     |  |  |  |  |  |
| 15. | 0°C Mini Cooler with Non toxic gel 32 places, Pack Size: 1                                                                                                                                                                                                            | 2  |     |  |  |  |  |  |
| 16. | 0°C Mini Cooler with Non toxic gel 12 places, Pack Size: 1                                                                                                                                                                                                            | 2  |     |  |  |  |  |  |
| 17. | PCR Mini Cooler, with Insulating solution 96 places, Pack Size: 1                                                                                                                                                                                                     | 2  |     |  |  |  |  |  |
| 18. | Ice Tray without Lid, 1L, Pack Size: 2                                                                                                                                                                                                                                | 2  |     |  |  |  |  |  |

|     |                                                                                                            |    |  |  |  |  |  |  |
|-----|------------------------------------------------------------------------------------------------------------|----|--|--|--|--|--|--|
| 19. | Ice Tray without Lid, 4L,<br>Pack Size: 1                                                                  | 2  |  |  |  |  |  |  |
| 20. | Parafilm 4"X 125', Pack Size: 1                                                                            | 1  |  |  |  |  |  |  |
| 21. | Petri dishes with triple vent<br>Radiation Sterile 90 X 14 mm,<br>Pack Size: 480 packed in sleeve<br>of 20 | 1  |  |  |  |  |  |  |
| 22. | Maxiamp PCR 330 Tubes 0.2<br>ml, Pack Size: 1000                                                           | 1  |  |  |  |  |  |  |
| 23. | Maxiamp 0.2 ml Tube Strips<br>with Cap 0.2 ml, Pack Size: 125                                              | 1  |  |  |  |  |  |  |
| 24. | PCR Workstation Pack 96<br>places, Pack Size: 6                                                            | 1  |  |  |  |  |  |  |
| 25. | Measuring Cylinder, 10<br>ml, Pack Size: 12                                                                | 1  |  |  |  |  |  |  |
| 26. | Measuring Cylinder, 100 ml,<br>Pack Size: 6                                                                | 1  |  |  |  |  |  |  |
| 27. | Measuring Cylinder, 1000 ml,<br>Pack Size: 4                                                               | 1  |  |  |  |  |  |  |
| 28. | Beaker with Handle, 1000 ML<br>Pack Size: 4                                                                | 1  |  |  |  |  |  |  |
| 29. | Beaker, 1000 ML<br>Pack Size: 4                                                                            | 1  |  |  |  |  |  |  |
| 30. | Analytical/Long Stem Funnel,<br>100 MM, Pack Size: 12                                                      | 1  |  |  |  |  |  |  |
| 31. | Racked Filter Tips Sterile 10 $\mu$ l,<br>Pack Size: 960 (10 racks)                                        | 2  |  |  |  |  |  |  |
| 32. | Racked Filter Tips Sterile 100<br>$\mu$ l, Pack Size: 960 (10 racks)                                       | 2  |  |  |  |  |  |  |
| 33. | Racked Filter Tips Sterile 200<br>$\mu$ l, Pack Size: 960 (10 racks)                                       | 2  |  |  |  |  |  |  |
| 34. | Racked Filter Tips Sterile 1000<br>$\mu$ l, Pack Size: 960 (10 racks)                                      | 2  |  |  |  |  |  |  |
| 35. | Microcentrifuge Tube, Sterile<br>1.5 ml, Pack Size: 100                                                    | 10 |  |  |  |  |  |  |
| 36. | Microcentrifuge Tube, Sterile<br>2.0 ml, Pack Size: 100                                                    | 10 |  |  |  |  |  |  |
| 37. | Float Rack, 16 places,<br>Pack Size: 6                                                                     | 1  |  |  |  |  |  |  |
| 38. | Float Rack, 8 places,<br>Pack Size: 6                                                                      | 1  |  |  |  |  |  |  |
| 39. | 96 Place Reversible Rack with<br>Cover, 96 places, Pack Size: 4                                            | 1  |  |  |  |  |  |  |
| 40. | Rack for Micro Tube, 20 places,<br>Pack Size: 6                                                            | 1  |  |  |  |  |  |  |
| 41. | Purepack Pasteur Pipette- sterile<br>3ml, Pack Size:                                                       | 1  |  |  |  |  |  |  |

|     |                                                                                                                    |   |  |  |  |  |  |  |
|-----|--------------------------------------------------------------------------------------------------------------------|---|--|--|--|--|--|--|
| 42. | Sterilized disposable flexible loop, 4.00mm diameter, calibrated to 0.01 ml individually sealed Pack Size: 100 No. | 1 |  |  |  |  |  |  |
| 43. | Pipettor Stand, capacity 5, Pack Size: 1 No.                                                                       | 2 |  |  |  |  |  |  |

### **Terms & Conditions.**

1. The interested agencies are required to submit the sealed technical bid and sealed financial bid separately. The bids in sealed Cover -I must comprise of "Technical Bid" and sealed Cover-II must comprise of "Financial Bid" Both these two cover must be placed in third sealed cover super-scribed with all the details mentioned (Quotation No/Name and Due date of opening must be written on top of the envelop) .
2. Technical Bid must comprise of complete specification of item and other related document as asked in the NIQ.
3. Sample will be called from the bidders who are technically qualified, if required.
4. Financial bids will be opened only to those bidders whose samples are technically qualified.
5. Rate should be mentioned in words & figures both.
6. GST rates applicable on the quoted item must be mentioned.
7. Delivery Schedule- within 15 days from the date of issue of P.O.
8. No additional documents related to this NIQ will be entertained after NIQ opened.
9. Price should be F.O.R. for Destination basis (i.e. Concerned Department).
10. LD@0.5% of delayed supply per week or part week for delay of supply of material subject to maximum upto 10% of delayed supply should be deducted.
11. The bidder participating in the LPC must have a local registered office in Chhattisgarh.
12. The GST registration details may please be furnished.

  
 4/10/25  
**PI - IRDL Project**  
**Department Microbiology**  
**AIIMS Raipur**  
**डॉ. अर्चना योगेन्द्र केचे**  
**Dr. Archana Yogendra Keche**  
 प्राध्यापक (सूक्ष्मजीव विज्ञान विभाग)  
 Professor (Department of Microbiology)  
 अखिल भारतीय आयुर्विज्ञान संस्थान, रायपुर (छ.ग.)  
 All India Institute of Medical Sciences, Raipur (C.G.)