



Date: - 19 /07/2021

Corrigendum

विषय/Sub: Corrigendum in Notice Inviting Tender for “Supply of Advanced Integrated Video Urodynamic Workstation” for Urology Department at AIIMS, Raipur

संदर्भ/Ref: Tender No.AIIMS/R/CS/URO/2020/32/OT

In context to following corrigendum is issued against subject tender:

Sr. No	Page/Clause/ Point	Existing Sentence	To be read as/ Addendum
01	Page 23, Point No.1 (1.11)	1.11. Nine Channel Esophageal Manometry	1.11. Nine Channel Esophageal Manometry (Optional)
02	Page 23, Point No.1 (1.12)	1.12. Near Infra-Red Spectroscopy (NIRS)	1.12. Deleted
03	Page 24, Point No.6 (6.2)	6.2. Should have facility to attach up to two Radio Frequency (RF) based, weight based Flow Transducer.	6.2. Should have facility to attach up to two Radio Frequency (RF)/ Bluetooth based, weight based Flow Transducer.
04	Page 24, Point No.6 (6.3)	6.3. Should be supplied with one wireless weight-based RF Uroflow Transducer and one wired weight-based Uroflow Transducer with flow range of 0-50 ml/sec.	6.3. Should be supplied with one wireless weight-based RF/ Bluetooth Uroflow Transducer and one wired weight-based Uroflow Transducer with flow range of 0-50 ml/sec.
05	Page 24, Point No.6 (6.4)	6.4. Flow transducer should be supplied with transparent security cover to avoid water spillage inside the transducer.	6.4. If required , Flow transducer should be supplied with transparent security cover to avoid water spillage inside the transducer.
06	Page 24, Point No.6 (6.8)	6.8. Should be supplied with Height-adjustable Micturation Chair with folding arms, detachable seat & cover - 1No., Removal Funnel - 2 Nos., Beaker/ Collection Jug – 2Nos., Transparent Security Cover for uroflowmetry Transducer -10 Nos.	6.8. Should be supplied with Height-adjustable Micturation Chair with folding arms, detachable seat & cover - 1No., Removal Funnel - 2 Nos., Beaker/ Collection Jug – 2Nos., If required , Transparent Security Cover for uroflowmetry Transducer - 10 Nos.
07	Page 24, Point No.9 (9.4)	9.4. Should be supplied with anal & vaginal probes.	9.4. Should be supplied with anal & vaginal probes/ Catheters for doing Biofeedback study.

08	Page 25, Point No.12	12. Should have Leak Point Detector for cough/ valsalva leak point pressure measurement.	12. Should have Leak Point Detection facility for cough/ valsalva leak point pressure measurement.
09	Page 25, Point No.14 (14.1)	14.1 Should be supplied with integrated mobile cart to fit the complete UDS along with All-in -One PC.	14.1 Should be supplied with integrated mobile cart to fit the complete UDS along with All-in -One PC/ Computer.
10	Page 25, Point No.17	17. Near Infrared Spectroscopy (NIRS): complete set - 1No.	17. Deleted
11	Page 26, Point No.19 (19.3)	19.3 Should have Display screen size 8-inch or more.	19.3 Should have Display screen size 6- inch or more.
12	Page 26, Point No.20 (20.1)	20.1 Double Lumen Cystometry Catheter 7-8Fr (Water filled/ Air Charged/ Microtip) - 10 Nos.	20.1 Double Lumen Cystometry Catheter 7-8Fr (Water filled/ Air Charged/ Microtip) - 30 Nos.
13	Page 26, Point No.20 (20.2)	20.2 Double Lumen Cystometry Catheter 5-6Fr (Water filled/ Air Charged/ Microtip) - 30 Nos.	20.2 Double Lumen Cystometry Catheter 5-6Fr (Water filled/ Air Charged/ Microtip) - 10 Nos.
14	Page 26, Point No.20 (20.13)	20.13 Transparent Security Cover for Uroflowmetry transducer - 10 nos.	20.13 If required , Transparent Security Cover for Uroflowmetry transducer - 10 nos.
15	Page 26, Point No.20 (20.14)	20.14 Anal Probe for biofeedback – 5 Nos.	20.14 Anal Probe/ Catheter for biofeedback – 5 Nos.
16	Page 26, Point No.20 (20.15)	20.15 Vaginal Probe for biofeedback – 5 Nos.	20.15 Vaginal Probe/ Catheter for biofeedback – 5 Nos.
17	Page 27, Point No.21 (21.10)	21.10 Should have US FDA and European CE with 4 digit notified body no or BIS approved certificate for the quoted model.	21.10 Should have US FDA or European CE with 4 digit notified body no or BIS approved certificate for the quoted model.

Note: All other terms & condition will be remain unchanged.

Stores Officer
AIIMS, Raipur (C.G.)